

---

**Particle Analysis and Display  
System (PADS):  
Cloud and Aerosol  
Spectrometer with  
Depolarization (CAS-DPOL)  
Module Manual**

**DOC-0275, Rev A**

**PADS 3.6.1, CAS-DPOL Module 3.6.1**



2545 Central Avenue  
Boulder, CO 80301 USA

*Copyright © 2011 Droplet Measurement Technologies, Inc.*

---

**2545 CENTRAL AVENUE  
BOULDER, COLORADO, USA 80301-5727  
TEL: +1 (303) 440-5576  
FAX: +1 (303) 440-1965  
WWW.DROPLETMEASUREMENT.COM**

All rights reserved. DMT licenses PADS software only upon the condition that you accept all of the terms contained in this license agreement. Each PADS license you purchase allows you to acquire data on one computer only. Data can be viewed in playback mode on an unlimited number of computers.

This software is provided by DMT “as is” and any express or implied warranties, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose are disclaimed. Under no circumstances and under no legal theory, whether in tort, contract, or otherwise, shall DMT or its developers be liable for any direct, indirect, incidental, special, exemplary, or consequential damages (including damages for work stoppage; computer failure or malfunction; loss of goodwill; loss of use, data or profits; or for any and all other damages and losses).

Some states do not allow the limitation or exclusion of implied warranties and you may be entitled to additional rights in those states.

### **Trademark Information**

All Droplet Measurement Technologies, Inc. product names and the Droplet Measurement Technologies, Inc. logo are trademarks of Droplet Measurement Technologies, Inc.

All other brands and product names are trademarks of their respective owners.

### **Risks of Installing Additional Software**

Instrument computers from DMT are configured to acquire data in a reliable, robust manner. Typically, such instruments are either not connected to a network or are connected to a small, local network that is isolated from the internet, reducing the risk of viruses. Since anti-virus programs can cause erratic behavior when run in the background on data acquisition computers, DMT does not install anti-virus, anti-spam, or anti-malware programs. If you choose to install these programs, you accept the risk associated with them in terms of potential performance degradation of the software installed by DMT.

For similar reasons, DMT recommends that you do not install or run other software on the dedicated instrument computer. Although the installation of some software may be unavoidable, it is particularly important not to run other software while the computer is acquiring data.

# CONTENTS

<b>1.0</b>	<b>Introduction .....</b>	<b>4</b>
<b>2.0</b>	<b>Configuration.....</b>	<b>4</b>
2.1	Configuring the CAS-DPOL.....	4
2.1.1	<i>CAS-DPOL Parameters</i> .....	5
2.2	Channels Table.....	7
2.3	Tables Tab .....	9
2.4	Configuring the CAS-DPOL Display .....	10
<b>3.0</b>	<b>The CAS-DPOL Window .....</b>	<b>11</b>
3.1	Sub-tabs.....	11
3.1.1	<i>Data Tab</i> .....	11
3.1.2	<i>LWC / # Conc and LWC / MVD Tabs</i> .....	12
3.1.3	<i>Selectable Charts</i> .....	12
3.1.4	<i>PBP Data</i> .....	12
3.1.5	<i>Tools Tab</i> .....	13
3.2	Selectable Channels Display .....	13
3.3	Histogram Data Window .....	13
<b>4.0</b>	<b>Zooming In and Out .....</b>	<b>14</b>
<b>Appendix A: CAS-DPOL Channels.....</b>		<b>15</b>
Main CAS-DPOL file .....		16
The CAS-DPOL PBP File .....		17
<b>Appendix B: Adjusting PBP Size Count Channels so They Scale Linearly .....</b>		<b>18</b>
<b>Appendix C: Revisions to Manual .....</b>		<b>19</b>

## ***List of Figures***

Figure 1:	CAS-DPOL Configuration Editor Window .....	5
Figure 2:	Example Channel Specifications in the Config Editor Window .....	8
Figure 3:	Tables Tab in CAS-DPOL Config Editor .....	9
Figure 4:	CAS-DPOL Display Editor Window.....	10
Figure 5:	Time-Range Controls .....	15

## ***List of Figures***

Table 1:	IPT Time Ranges for CAS-DPOL IPT Bins .....	12
----------	---	----

## 1.0 Introduction

The Particle Analysis and Display System (PADS) is a software package that interfaces with instruments produced by Droplet Measurement Technologies (DMT) and other leading instruments used in the atmospheric sciences. This manual describes the PADS Cloud and Aerosol Spectrometer with Depolarization (CAS-DPOL) module version 3.6.1.

For an explanation of the basic PADS setup and instructions on how to acquire data using PADS, consult the *PADS Overview Manual, DOC-0300*. Definitions and calculations used in the CAS-DPOL module are also described in the *PADS Overview Manual*.

## 2.0 Configuration

Using PADS, you can configure both the software settings for the instrument and the instrument's data display in PADS. The following two sections explain how to do this. Configuring the instrument's software and display affects the default settings PADS uses when starting up. Some parameters can also be changed while PADS is running, but doing so affects the current session only.

### 2.1 Configuring the CAS-DPOL

Your CAS-DPOL and data system should arrive preconfigured from DMT. In some cases, *however, you may want to change the software configuration for the instrument. To do this, follow the steps below. Note: Droplet Measurement Technologies STRONGLY recommends that customers contact our office prior to changing any of the parameters in the instrument configuration. Improper changes can result in communication failure and/or changes in PADS computation algorithms, which can compromise data validity.*

1. Click on the "CAS-DPOL" tab.
2. From the **Configure** menu, select **Configure Instrument**. You will see the following window.




Figure 1: CAS-DPOL Configuration Editor Window

3. Now you can configure the instrument parameters to your desired specifications. See the definitions below for explanations of individual parameters. If at any time you would like to revert to the previously saved values for the CAS-DPOL parameters, press **Cancel** to exit the window without saving changes. Pressing **Reset Parameters** reverts parameters to their DMT-supplied default values.
4. When you are done configuring the CAS-DPOL parameters, press **Save** at the top of the Config editor window. Then press the green **Reset Program** button for the new configuration to take effect. Note that pressing the **Reset Program** button will clear any data currently being displayed.

### 2.1.1 CAS-DPOL Parameters

**Instrument #:** This lists the number corresponding to the instrument you are viewing, in this case the CAS-DPOL. If your CAS-DPOL has been assigned instrument number one, you will see “1” in this field. You should not need to modify the instrument number, and in fact you are unable to do so from within PADS.

**Sample Time:** This parameter shows the time interval you’d like between samples. You can have the probe sample at intervals between 0.04 and 20 sec (25 to 0.05 Hz). Note that if you increase the sample time, you will still collect data for the same number of particles. This is because the probe collects data continuously and relays cumulative data at each sampling

interval. For example, say you have the sample time set to .5 seconds. You might see four particles of size 25  $\mu\text{m}$  during the first sample, and five particles of this size during the second sample. If you had set your sample time to one second instead of .5 seconds, you would instead get one sample that showed nine particles of size 25  $\mu\text{m}$ . *Note:* Sample Time is most often set to 1 Sec. Higher sample rates may or may not work on a given data system, depending on the computer performance and the number and types of instruments PADS is configured to use.

**At Startup Enabled / Disabled:** If you want the CAS-DPOL to acquire data when PADS begins sampling, make sure this parameter is in the “Enabled” mode. In some cases, such as if the CAS-DPOL is inoperative, you may want to use this control to disable the probe. Disabling the CAS-DPOL allows data to transmit from other instruments without interference. Data will still be written to the disabled instrument’s output file, but PADS will write “NaN” to all fields.

**COM Port:** This is the serial communications port that the CAS-DPOL uses to connect with the computer. This number should match the computer hardware configuration for the particular computer you are using. If you are not using multiple computers, this number should not be changed.

**Baud Rate:** The baud rate for the probe is defined at manufacture. This parameter has been grayed out and you should not need to change it. If you reconfigure your hardware, however, the baud rate may change. If this occurs, contact DMT for help in changing your baud rate in PADS.

**FW Rev.:** This parameter lists the firmware revision number. This parameter has been configured to work for your particular instrument and should not be changed. PADS uses the firmware revision number when calculating the **Back : Fwd** and **Dpol : Back** ratios. See Appendix B for details.

**Sample Area:** This is the physical area in which particles are detected. CAS-DPOL sample area is a constant regardless of particle size. This value is preconfigured to match your instrument, so it is strongly suggested that you do not change it.

**Channel Count:** This number indicates how many sizing bins the CAS-DPOL uses to categorize particles. This number has been grayed out because it is preconfigured for your instrument and should not change.

**Recovery Factor:** This parameter is used in calculating ambient temperature from measured temperature. PADS uses Bernoulli's equation for this calculation. For more information on this equation and the recovery factor, see the "Ambient Temperature" entry in Appendix B of the *PADS Overview Manual*. By default, Recovery Factor is set to 1.0.

The **PAS Source** control specifies from which of the following sources the system should obtain the applied probe air speed (PAS):

- 1.) A specific instrument in the system (this can be any instrument capable of measuring air speed)
- 2.) A manually entered value:
  - a. A "Local" value, which at start-up is the value entered in the box to the right of the source control. This number can be changed from the instrument display while the program is running.
  - b. A "Global" value entered on the **Setup** tab

Applied PAS is used to calculate sample volume. In flight conditions, you will typically want to select an instrument as the air speed source. However, you will need to enter manual air speed values during probe calibration.

Pressing the **Reset Parameters** button resets all parameters to their DMT-supplied default values. After making changes to the parameters, you will need to press the **Save** button and then click the green **Reset Program** to activate these changes. Clicking **Reset Program** will clear any data PADS is currently displaying.

## 2.2 Channels Table

The channels listed in the Channels table are configurable. These are A/D housekeeping channels that measure a 0 - 10 V range from one of the instrument's internal sensors, for example a pressure or temperature sensor. A conversion equation converts the A/D counts into other, more meaningful units (e.g., mBar or °C). You can specify this equation in the Channels table.

Note: While it is possible to use the Channels table to rename output channels, in most cases your system is preconfigured so that the channels in the table correctly match output from your instrument(s). While minor rescaling of output channels can improve data accuracy, DMT does not recommend altering your basic channel configuration.

The second column in the table indicates the type of equation that PADS should use—linear, polynomial, or none. (“Thermister D” and “Thermister G” are complicated, pre-set equations used by some instruments, while “Custom” allows users to select an equation they have entered on the Configure menu.) “Linear” indicates a linear equation, while “4<sup>th</sup> Order Poly” indicates a higher order polynomial equation with up to five terms. “None” means the digital value (between 0 and 4095) will be returned without any scaling.

The right-hand fields in the channels table indicate the coefficients to be used in the conversion equation. Figure 2 shows the setup for a hypothetical channel with the second-order polynomial conversion equation, as follows:

$$\text{New\_Channel} = 34.01 + 0.061 x + 0.0092 x^2$$

where *x* is the digitized analog value returned by the A/D converter.




Figure 2: Example Channel Specifications in the Config Editor Window

The number of coefficients that PADS uses depends on the equation type. “None” does not use any coefficients. “Linear” uses the first two coefficients, which are listed in the first two table cells after the equation type. “4th Order Poly” uses one to five coefficients.

In cases where there are non-zero numbers in cells that are not used in the function, PADS ignores these numbers. For instance, if you specify “Linear” as your function and have .32 in the farthest right cell, the program will just ignore this value.

Clicking the Reset Channels knob at the bottom of the CAS-DPOL Parameter window will reset all the channel parameters to their DMT-supplied default values.

After making changes to the Channels tab, you will need to press the Save button and then click the green Reset Program to activate these changes. Clicking Reset Program will clear any data PADS is currently displaying.

## 2.3 Tables Tab

The **Tables** tab lists the threshold tables that allow the CAS-DPOL to bin particles according to size.




Figure 3: Tables Tab in CAS-DPOL Config Editor

**Bin 1 Lower Thresh.** lists the lower size limit in  $\mu\text{m}$  of the smallest bin. **ADC Lower Thresh.** indicates the smallest peak A/D value a particle can have and still be sized by the instrument. Small noise peaks can occasionally occur in the absence of particles, and the ADC Lower Thresh. can be used to eliminate such noise. Both of these parameters can be changed by typing a new value into the text box or by using the control arrows to the left.

The table below consists of five columns. The first column lists the bin number. The second and third columns are for the forward-scattering optics. The second column lists the bin’s upper size limit in microns, while the third column lists the A/D count that corresponds to this upper size limit. Lower bin boundaries are the upper limit of the previous bin, except for bin 1, whose lower boundaries are listed above the table.

Bin numbers are fixed, but you can modify both the upper size boundaries and the corresponding A/D count values. To modify a cell, either type in a new value, or use the arrows to the left. The scrollbar to the right allows you to access other rows in the table.

If the threshold tables contain incorrect values, PADS will alert you by displaying a yellow error message. Errors will result if values are not in ascending order or if there is not a value for each bin number listed in column one. If there are extra rows in the table, you can remove them by clicking on the **Remove extra array indices** button. This button only appears if extra indices exist.

The **Reset Tables** knob restores table values to their DMT-supplied defaults.

After making changes to the Tables tab, you will need to press the **Save** button and then click the green **Reset Program** to activate these changes. Clicking **Reset Program** will clear any data PADS is currently displaying.

## 2.4 Configuring the CAS-DPOL Display

To configure the CAS-DPOL display, click on the CAS-DPOL tab if you have not already done so. Then select **Configure** from the menu bar and click on **Configure Display**. This will bring up the following window.




Figure 4: CAS-DPOL Display Editor Window

You do not need to modify the **Display #** or **Instrument #**.

Changing the **Refresh Time** allows you to set the time intervals for data display during acquisition mode; you can choose any time that is equal to or greater than the sample time. (Choosing a time less than the sample time is not useful, since the same data will be displayed multiple times.)

**Chart 1** and **Chart 2** allow you to configure the channels on the CAS-DPOL's selectable graphs. To change these channels, click on the arrow buttons to scroll between available options for the channels. You can also click on the white fields to bring up a list of all the available options, from which you can then choose the channel you want.

**Histogram Label:** This control allows you to specify the units for the x-axis of the histogram, which can either be size in  $\mu\text{m}$  or bin numbers.

The **Selectable Channels Display** controls which channels PADS displays in the upper right of the CAS-DPOL window.

When you are done, click on **Save** to update the configurations or **Cancel** to revert to the previous configuration. After you reset PADS, you will be able to see any changes. Note that clicking **Reset Program** will clear out any data currently being displayed.

Configuring channels in the **Display Editor** will change the display upon start-up. Once PADS has started, you can change many of these settings from within the main CAS-DPOL tab.

## 3.0 The CAS-DPOL Window

The different parts of the CAS-DPOL Window are discussed below. For explanations of the **Enable** button, **COM Port** indicator, and **Fault/No Fault** button, see the “Instrument Tabs” section of the *PADS Overview Manual*.

### 3.1 Sub-tabs

The CAS-DPOL window has six sub-tabs, which are discussed below.

#### 3.1.1 Data Tab

The **Data** tab displays the current values for many of the CAS-DPOL data channels. For more information on these channels, see Appendices A and B of the *PADS Overview Manual*.

### 3.1.2 LWC / # Conc and LWC / MVD Tabs

The **LWC / # Conc** tab displays a time-trace chart of calculated liquid water content in  $\text{g}/\text{m}^3$  (shown in green) and # Conc in particles/ $\text{cm}^3$  (shown in blue). The **LWC /# MVD** tab displays a time-trace chart of calculated liquid water content in  $\text{g}/\text{m}^3$  (green) and MVD in  $\mu\text{m}$  (blue). A red cursor shows the current moment in time.

### 3.1.3 Selectable Charts

This tab displays time-trace charts of two user-selectable channels. The charts are overlaid on each other, with the left axis (in green) corresponding to the channel listed in the above left of the chart, and the right axis (in blue) applying to the channel listed in the above right. To change these channels, click on the names of the currently displayed channels, which will bring up a list of options. If you would like to permanently change the **Selectable Charts** channels, as opposed to simply changing them for the current session, you can do so by selecting **Configure > Configure Display** and making the desired changes.

### 3.1.4 PBP Data

The PBP data tab shows particle-by-particle data. On the right are various ratios involving the backscattered signal, forward-scattered signal, and depolarization signal. These ratios are averages for the sampling instance. The tab also displays the average interparticle time (PBP Avg IPT) in msec and the average standard deviation of this time, also in msec.

The histogram shows particle counts binned by IPT time. IPT Time ranges for each bin are given in Table 1. Note that the Bin 28's upper boundary is really 1677.72, not the 1000 that PADS displays. This is the longest inter-particle time that the CAS-DPOL can detect.

Bin	IPT (msec)	Bin	IPT (msec)	Bin	IPT (msec)	Bin	IPT (msec)
1	0 - 1	8	7 - 8	15	50 - 60	22	300 - 400
2	1 - 2	9	8 - 9	16	60 - 70	23	400 - 500
3	2 - 3	10	9 - 10	17	70 - 80	24	500 - 600
4	3 - 4	11	10 - 20	18	80 - 90	25	600 - 700
5	4 - 5	12	20 - 30	19	90 - 100	26	700 - 800
6	5 - 6	13	30 - 40	20	100 - 200	27	800 - 900
7	6 - 7	14	40 - 50	21	200 - 300	28	900 - 1678

Table 1: IPT Time Ranges for CAS-DPOL IPT Bins

The **Y-scale** button allows you to scale histogram data linearly or logarithmically.

### 3.1.5 Tools Tab

The Tools tab lists various **Source** controls. If these controls are grayed out, click on **Press to Enable Source Changes** to enable them.

The **PAS Source** control specifies from which of the following sources the system should obtain the applied probe air speed (PAS):

- 1.) A specific instrument in the system (this can be any instrument capable of measuring air speed)
- 2.) A manually entered value:
  - a. A “Local” value, which at start-up is the value entered in the box to the right of the source control. This number can be changed from the instrument display while the program is running.
  - b. A “Global” value entered on the **Setup** tab

Applied PAS is used to calculate sample volume. In flight conditions, you will typically want to select an instrument as the air speed source. However, you will need to enter manual air speed values during probe calibration.

Changing the source settings here changes them for the current session only. (To change the default values used upon start-up, do so from the **Configure > Configure Instrument** menu.) After you have finished specifying new sources, click on the **Press to Lock Source Changes** button to disable further changes.

*Note:* the PAS source can be changed during acquisition, but if data are reprocessed, those changes to the source that occurred during acquisition will be lost; see the “Configure Menu” section of *DOC-0300*, the *PADS Overview Manual*.

## 3.2 Selectable Channels Display

The upper right of the CAS-DPOL window displays current data for ten user-selectable channels. New channels can be selected by clicking on the name of a current channel and selecting a replacement from the drop-down list. To make permanent changes to this list, edit the settings on the CAS-DPOL Display Editor (select **Configure > Configure Display**).

## 3.3 Histogram Data Window

Below the CAS-DPOL tabs is the histogram display of particle data. The histogram shows the particle distribution detected by the forward-scattering optics. The number in the upper right of the histogram shows the total particle count across the entire histogram.

The x-axis labels indicate the upper boundary of each bin. A bin's lower boundary is the upper boundary of the previous bin, except for the first bin. You can view and modify the lower boundary of this bin—and other bins—by selecting **Configure > Configure Instrument** and then clicking on the **Tables** tab.

To the right of the histogram are buttons that control the scaling and display of the histogram data. **Auto-scale** controls let you control how the x and y scales are set. If you enable autoscaling, PADS will automatically select an appropriate scale with which to display the current data. For instance, if the probe is not currently detecting many particles, the y-axis range will decrease. On the other hand, if you disable autoscaling, the scale of the axes will remain constant. In this case, the range will always be the same as it was when autoscaling was disabled. The minimum and maximum values can then be changed manually by typing new numbers into these fields.

When you turn on the **Log-Scale** buttons, PADS scales the appropriate axis logarithmically rather than linearly. Autoscaling can be enabled or disabled with this option.

If you click on the **Normalized** button, PADS will scale the particle data so that each bin of data is divided by the width of that bin in  $\mu\text{m}$ . Normalization is useful when the widths of the size bins are not constant, as is often the case when an optical spectrometer has to cover a large size range or if some size bins are made larger to decrease uncertainty due to operational limitations. Without normalization, the concentrations from one size bin cannot be compared quantitatively from those in another size bin of different width. For example, if 8 particles are counted in a size bin that is 2  $\mu\text{m}$  wide and 8 particles are counted in the next larger channel that is 4  $\mu\text{m}$  wide, the size distribution would appear flat and we would assume that it was equally probable to measure particles of either size. This is a biased sample, however, since the larger size bin is twice as wide as the smaller. If we normalize by the width, however, we find that the smaller size category has 4 particles per  $\mu\text{m}$  and the larger size bin only has 2 particles per  $\mu\text{m}$ . This is now an unbiased sample.

The **size / bins** radio buttons allow you to change the x-axis units. When the x-axis displays bin number, each bin's upper size limit in  $\mu\text{m}$  is displayed above each data bar.

## 4.0 Zooming In and Out

There are several ways to zoom in or out on CAS-DPOL charts and the histogram. As described in the *PADS Overview Manual*, you can use the time-range controls (Figure 5) to zoom. To zoom in on the data, move the green and red controls close to the white control, which will narrow the range of displayed data. To zoom out, move the two colored controls away from the white control.




Figure 5: Time-Range Controls

On the chart itself, you can also type numbers directly into the first and last labels on the x and y axis to change the scaling.

*Note:* Do not right-click on chart and change the auto-scaling using the drop-down menu. This can interfere with the chart display. PADS autoscales most charts automatically. You can turn off autoscaling on the histogram using the buttons in the lower right of the window.

## Appendix A: CAS-DPOL Channels

The CAS-DPOL generates two processed output files. These will be named something like the following:

- 00CAS DPOL20110311090831.csv
- 00CAS DPOL20110311090831\_PBP.csv

The first file contains CAS-DPOL output gathered over each sampling instance. This file contains data on all the particles the probe has detected. The second file, which has “\_PBP” appended to its name, contains particle-by-particle data. PADS collects particle-by-particle data only for the first 292 particles detected in each sampling instance. If more than 292 particles arrive in the sampling interval, they are ignored in the second file but included in calculations of concentration, LWC, and MVD that are recorded in the first file. If no qualified particles are detected during a given sampling instance, the PBP file will not store any data. More information about these two files is given below.

CAS-DPOL channels fall into several broad categories:

- Time channels
- Bin channels, which store data on the number of particles of different sizes that the CIP has detected
- Particle statistics (e.g., rejected particles, number concentration, etc.)
- Probe statistics
- Channels reserved for internal use (e.g., Host Sync Counter and Reset Flag)
- Particle-by-particle channels

## Main CAS-DPOL file

The main file will typically have the channels listed below. The channels in italics are configurable housekeeping channels, so they may be different on your instrument. Spare channels 1 - 7 are generally not used on CAS-DPOLs that are part of multi-instrument systems. On stand-alone CAS-DPOLs, these channels are renamed and store data relating to pressure, temperature, and the photodiodes (see the list below) You can plot any CAS-DPOL channel with respect to time using the Selectable Charts tab.

A list of channels for a 30-bin CAS-DPOL appears on the following page. The CAS-DPOL output file will contain data values for each channel for each sampling instance. You can also plot each of these channels with respect to time using the CAS-DPOL Selectable Charts tab. For definitions of the channels, consult *Appendix A: Definitions* in the *PADS Overview Manual*.

End Seconds	<i>Back High Gain Baseline (V)</i>
Day of Year	<i>Spare</i>
Year	<i>Dpol Baseline (V)</i>
Status	<i>Electronics Temp (C)</i>
Sum of Transit	<i>Dpol TEC Temp (C)</i>
Sum of DOF	<i>Dpol Heat Sink Temp (C)</i>
Fwd OverRange	<i>LWC Hotwire (V)</i>
Back OverRange	<i>LWC Slave (V)</i>
DPOL OverRange	<i>Laser Current (mA)</i>
Spare 1 - 3	<i>Laser Monitor</i>
<i>Dynamic Press (mbar)</i>	Spare 1 - 8
<i>Static Press (mbar)</i>	T Ambient (C)
<i>Recovery Temp (C)</i>	PAS (m/s)
<i>Fwd Heat Sink Temp (C)</i>	Number Conc (#/cm <sup>3</sup> )
<i>Back Heat Sink Temp (C)</i>	LWC (g/cm <sup>3</sup> )
<i>Fwd Block Temp (C)</i>	MVD (um)
<i>Back Block Temp (C)</i>	ED (um)
<i>Spare 4 - 7</i>	Back : Fwd
<i>Qual TEC Temp (C)</i>	Dpol : Fwd
<i>Fwd TEC Temp (C)</i>	Dpol : Back
<i>Back TEC Temp (C)</i>	PBP Avg IPT (msec)
<i>Qual Heat Sink Temp (C)</i>	PBP IPT Std Dev (msec)
<i>Qual High Gain Baseline (V)</i>	Applied PAS (m/s)
<i>Qual Mid Gain Baseline (V)</i>	CAS Bin 1 – 30
<i>Qual Low Gain Baseline (V)</i>	CAS IPT Bin 1 – 28
<i>Fwd High Gain Baseline (V)</i>	UTC Seconds / GPS Time
<i>Fwd Mid Gain Baseline (V)</i>	Date
<i>Fwd Low Gain Baseline (V)</i>	Time

If there is no instrument in the system that reports **GPS Time**, or if such an instrument exists but the user has selected on the CAS-DPOL Config Editor to show UTC Seconds, the output channel file will contain **UTC Seconds**. Otherwise, it will report **GPS Time**.

The last two channels, **Date** and **Time**, will be listed after **GPS Time/ UTC Seconds** only if **Write Date & Time Stamp** is enabled on the **Setup** tab.

Many of the probe statistics are stored in “housekeeping channels,” a term that refers to data gathered with A/D sensors. The CAS-DPOL has 31 A/D housekeeping channels that have a 0-10 V range measured by a 12-bit A/D converter that gives integer values from 0 to 4095. Several of the housekeeping channels store data that indicate whether the probe is functioning properly. Several others are non-functional. Housekeeping channels are denoted by italics in the list above.

## The CAS-DPOL PBP File

The output second file has **\_PBP** appended to its name and contains particle-by-particle data. The file lists the following channels:

End Seconds  
Forward Size [counts]  
Back Size [counts]  
DePol Size [counts]  
IPT [msec]  
Date  
Time

The last two channels, **Date** and **Time**, are listed only if **Write Date & Time Stamp** is enabled on the **Setup** tab.

These channels get reported for the first 292 particles that the CAS-DPOL detects in each sampling instance.

## Appendix B: Adjusting PBP Size Count Channels so They Scale Linearly

The CAS-DPOL pbp output file has three channels that store data on particles' peak signals: **Forward Size**, **Back Size**, and **DePol Size**. Forward Size does not scale linearly, nor does Back Size on older versions of CAS-DPOL firmware (version 1.0.x). This non-linear scaling occurs because the forward-scattering detector and backward-scattering detector can use different gain stages for particles that emit weaker or stronger signals. On firmware version 2.0.x, Back Size does scale linearly. DePol Size always scales linearly, regardless of firmware version. (Note that linearity in this discussion refers only to optical-scattering cross-section measurements, not to particle size.)

In order to calculate the ratios **Back : Fwd**, **Dpol : Fwd**, and **Dpol : Back**—three channels that appear in the main CAS-DPOL output file—PADS must first scale Forward Size and Back Size (if applicable) linearly. It uses the algorithms given below to do so. These algorithms are also useful should you want to convert Forward and Backward PBP data to a linear scale for your own analysis.

Note that there is a small overlap between gain stages, which depends on the calibration of your particular instrument. These calibrations are documented for your instrument, though accurate interpretation of them may require additional consultation with DMT. As a result of the overlapping gain stages, even the adjusted scales will not be perfectly linear; where overlap occurs, digital counts may jump slightly. Nonetheless, the adjusted sizes yielded by the algorithms correspond much more closely to scattered light than do the raw forward and backward size channels.

### *To Scale Forward Sizes - For All Firmware Versions:*

Assume [Forward Size] is your raw reading and you wish to obtain an adjusted, linearly scaled reading, [Adjusted Forward]:

If  $20 \leq [\text{Forward Size}] \leq 3071$ ,

$$[\text{Adjusted Forward}] = [\text{Forward Size}]$$

If  $3072 \leq [\text{Forward Size}] \leq 6143$

$$[\text{Adjusted Forward}] = ([\text{Forward Size}] - 3071) \bullet 22$$

If  $6143 < [\text{Forward Size}] < 9216$

$$[\text{Adjusted Forward}] = ([\text{Forward Size}] - 6143) \bullet 506$$

The values 22 and 506 represent the multiplication factors between the mid and low gain stages, respectively, and the high gain stage. If [Forward Size] exceeds 9216, the particle is oversized and PADS will not size it in the histogram. Instead, it will be counted in the Fwd\_Overflow channel in the main CAS-DPOL output file. If the particle has been observed during particle-by-particle analysis, the Forward Size channel in the PBP file will be set to 12287.

#### *To Scale Backward Sizes - For Firmware Version 1.0.x:*<sup>1</sup>

Assume [Backward Size] is your raw reading and you wish to obtain an adjusted, linearly scaled reading, [Adjusted Backward]:

If  $0 \leq [\text{Backward Size}] \leq 1536$ ,

$$[\text{Adjusted Backward}] = [\text{Backward Size}]$$

If  $1537 \leq [\text{Backward Size}] \leq 3071$

$$[\text{Adjusted Backward}] = ([\text{Backward Size}] - 1536) \bullet 22$$

If [Backward Size] exceeds 3071, PADS counts it in the Fwd\_Overflow channel in the main CAS-DPOL output file. If the particle has been observed during particle-by-particle analysis, the Back Size channel in the PBP file will be set to 4095.

Note: the calibration coefficients 22 and 506 are the ones PADS uses to compute the three ratio channels mentioned above. As discussed, however, the coefficients for your particular instrument may be slightly different. For more information, consult DMT.

## **Appendix C: Revisions to Manual**

This manual replaces DOC-0189, the *CAS-POL-PBP PADS Operator Manual* for PADS version 2.X. All sections have been updated.

---

<sup>1</sup> Firmware version 2.0.x does not require scaling of the back-scatter signal.