




PARTICLE ANALYSIS AND DISPLAY SYSTEM (PADS)



SPECIFICATIONS:

- Flexible software for data acquisition and playback
- Interfaces with all DMT instruments
- Easy configuration of new instruments
- Runs under Windows 2000 or Windows XP
- Data files can be imported into data base or graphing programs for additional analysis
- 0.1 to 10 Hz sampling speed
- Tab structure for rapid view of multiple instruments
- Unlimited distribution of playback program
- Separate program for review of 2D images
- Can be supplied as the software only or with a completely configured data system